

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cinnamon leaf oil - Ceylon origin, FG

Product Number: **W229210**

CAS Number: 8015-91-6

TEST	Specification
Appearance (Color)	Yellow to Brown
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.529 - 1.537
Infrared spectrum	Conforms to Structure
Specific Rotation	-2 - 1 °
Assay	75.00 - 88.00 %
Specific Gravity	1.030 - 1.050
Solubility (Turbidity)	Clear
Solubility (Color)	Yellow to Brown
1ml/1.5ml 70% Ethanol	
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.2.ZQ5.10000029935

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.