


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Gallic acid - 97.5-102.5% (titration)

Product Number: G7384
CAS Number: 149-91-7
MDL: MFCD00002510
Formula: C₇H₆O₅
Formula Weight: 170.12 g/mol



TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Very Slightly Hazy
50 mg/mL, EtOH	
Loss on Drying	≤ 10 %
Titration by NaOH (dry basis)	97.5 - 102.5 %
Purity (GC)	≥ 98.5 %

Specification: PRD.1.ZQ5.10000011445

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.