

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Dodecylbenzenesulfonic acid solution - 70 wt. % in isopropanol

Product Number:

522953

CAS Number:

27176-87-0

MDL:

MFCD00066677

Formula:


C₁₈H₃₀O₃S

Formula Weight:

326.49 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

70.0 - 80.0 %

Specification: PRD.0.ZQ5.10000055819

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.