


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Methoxy-4-nitrobenzoic acid - 98%

Product Number: 184306
CAS Number: 5081-36-7
MDL: MFCD00007353
Formula: C₈H₇NO₅
Formula Weight: 197.14 g/mol



TEST	Specification
Appearance (Color) Pale Yellow to Yellow or Tan	Conforms to Requirements
Appearance (Form) Powder, Crystalline Powder or Crystalline Fibers	Conforms to Requirements
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 97.5 %

Specification: PRD.0.ZQ5.10000042892

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.