

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

1-(4-METHOXYPHENYL)PIPERAZINE, 97%

Product Number:

571415

CAS Number:

38212-30-5

MDL:

MFCD00040736

Formula:

C₁₁H₁₆N₂O

Formula Weight:

192.26 g/mol

TEST

Specification

Appearance (Color)

White to Faint Brown and Faint Beige

Appearance (Form)

Powder or Crystal or Chunk(s)

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000039563

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.