

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Phosphorus tribromide - 97%

Product Number: 157783
CAS Number: 7789-60-8
MDL: MFCD00011436
Formula: Br₃P
Formula Weight: 270.69 g/mol



TEST	Specification
------	---------------

Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless to Yellow and Colorless to to Light Brown-Yellow
Appearance (Form)	Liquid
Titration by AgNO ₃	96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000069704

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.