

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Platinum on alumina - extent of labeling: 1 wt. % loading, powder

Product Number:

205966

Pt

MDL:

MFCD00011179

TEST	Specification
------	---------------

Appearance (Color)

Conforms to Requirements

Brown to Grey

Appearance (Form)

Powder

Moisture

Specification: Report Results

ICP Major Analysis

Confirmed

Confirms Aluminum Component

Assay

0.70 - 1.30 %

Platinum by ICP

Specification: PRD.2.ZQ5.10000032016

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.