

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
5,6,7,8-Tetrahydroquinoxaline - ≥97%, FG

Product Number: W332100
CAS Number: 34413-35-9
MDL: MFCD00038078
Formula: C₈H₁₀N₂
Formula Weight: 134.18 g/mol



TEST

Specification

| | |
|----------------------------|-------------------------------|
| Appearance (Color) | Colorless to Yellow to Orange |
| Appearance (Form) | Liquid |
| Refractive index at 20 ° C | 1.540 - 1.547 |
| Infrared spectrum | Conforms to Structure |
| Purity (GC) | ≥ 97.0 % |
| Impurities | ≤ 0.09 % |
| Aniline | |
| Arsenic (As) | ≤ 3 ppm |
| Cadmium (Cd) | ≤ 1 ppm |
| Mercury (Hg) | ≤ 1 ppm |
| Lead (Pb) | ≤ 10 ppm |
| Expiration Date Period | ----- |
| 5 Years | |

Specification: PRD.1.ZQ5.10000064392

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.