


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-(Trifluoromethyl)phenylacetic acid - 97%

Product Number: 233021
CAS Number: 32857-62-8
MDL: MFCD00004352
Formula: C₉H₇F₃O₂
Formula Weight: 204.15 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Pale Yellow
Appearance (Form)	Powder or Crystals or Chunk(s)
Infrared spectrum	Conforms to Structure
Proton NMR spectrum	Conforms to Structure
Purity (Titration by NaOH)	96.5 - 103.5 %
GC (area %) (VPCT)	≥ 96.5 %

Specification: PRD.0.ZQ5.10000064904

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.