

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Acetal - ≥98%, FG

Product Number: W200204
CAS Number: 105-57-7
MDL: MFCD00009243
Formula: C₆H₁₄O₂
Formula Weight: 118.17 g/mol



TEST

Specification

Appearance (Color)	Colorless
Appearance (Form)	Liquid
Refractive index at 20 ° C	1.379 - 1.383
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.0 %
Color Test	≤ 100 APHA
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Trace Metal Analysis	≤ 10 ppm
Expiration Date Period	-----
5 Years	

Specification: PRD.0.ZQ5.10000019860

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.