

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Hemoglobin A<sub>2</sub>, Ferrous Stabilized human - lyophilized powder

Product Number: **H0266**  
CAS Number: 9034-53-1  
MDL: MFCD00131286

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	Red
Appearance (Form)	Powder
Solubility (Color)	Red
Solubility (Turbidity)	Clear
20 mg/mL, H <sub>2</sub> O	
% Protein Drabkin	40 - 55
Approx. 50%	
Methemoglobin	≤ 20 %
Agarose Electrophoresis	97 - 100 %
Approx. 98%	
Tested For Infectious Agents	Tested

Specification: PRD.0.ZQ5.10000049801

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.