

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Nitric Oxide Synthase, Endothelial (1185-1205) antibody produced in rabbit - IgG fraction of antiserum, buffered aqueous solution

Product Number:

N3893

TEST**Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Colorless
IgG Concentration (OD 280nm)	9.0 - 15.0 mg/ml
WB-Cell Line/Tissue Extract	Bovine Lung Endothelial Cells
WB Titer	1:10,000-1:20,000
IHC - Tissue	Mouse Heart
IHC - Fixation	FFPE+ Antigen Retrieval-Decloaking or Frozen Sections Acetone Fix
IHC - Titer	1:100-1:200
Recommended Retest Period	-----
4 Years	

Specification: PRE.0.ZQ5.10000020528

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.