

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Samarium(II) iodide solution - 0.1 M in THF, contains samarium chips as stabilizer

Product Number:

347116

MDL:

MFCD00058873

Formula:

I₂Sm

Formula Weight:

404.17 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Dark Blue to Very Dark Blue

Appearance (Form)

Liquid

Purity Titration with Iodine

≥ 0.075 M

Expiration Date Period

9 Months

Specification: PRD.5.ZQ5.10000004226

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.