

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium 2-fluoro-5-formylphenyltrifluoroborate

Product Number:

723991

CAS Number:


1012868-70-0

Formula:

C₇H₄BF₄KO

Formula Weight:

230.01 g/mol

**TEST****Specification**

Appearance (Color)

White to Beige

Appearance (Form)

Powder or Crystals

Proton NMR spectrum

Conforms to Structure

Fluorine NMR Spectra

Conforms to Structure

Carbon

32.1 - 41.0 %

Purity (HPLC)

≥ 88.0 %

ICP Major Analysis

Confirmed

Confirms B and K Components

Specification: PRD.1.ZQ5.10000005671

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.