

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
DL-Panthenol - 99%

Product Number: 295787
CAS Number: 16485-10-2
MDL: MFCD00002944
Formula: C₉H₁₉NO₄
Formula Weight: 205.25 g/mol



TEST

Specification

Appearance (Color)

White

Appearance (Form)

Powder or Chunks

Infrared spectrum

Conforms to Structure

GC (area %)

≥ 98.5 %

Recommended Retest Period

6 Years

Specification: PRD.1.ZQ5.10000010255

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.