


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ammonium sulfate – suitable for plant cell culture, ≥99.0%

Product Number: A3920
CAS Number: 7783-20-2
MDL: MFCD00003391
Formula: H₈N₂O₄S
Formula Weight: 132.14 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
200 g + 300 ml, H ₂ O	
Heavy Metals <= 2ppm (as Lead)	Pass
Purity (Titration by NaOH)	≥ 99.0 %
Recommended Retest Period	-----
5 years	
Plant Cell Culture Test	Pass

Specification: PRD.1.ZQ5.10000047800

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

