


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
1-Dodecanol – Vetec™ reagent grade

Product Number: V900234
CAS Number: 112-53-8



Formula: C₁₂H₂₆O
Formula Weight: 186.34 g/mol

TEST

Specification

Appearance (Colour)	Conforms to Requirements
Appearance (Form)	Conforms to Requirements
Solid or Crystalline Solid or Liquid	
Purity (GC)	≥ 97.5 %
Infrared spectrum	Conforms to Structure

Specification: PRE.0.ZQ5.10000002515

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.