

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Poly(trimellitic anhydride chloride-co-4,4'-methylenedianiline) - powder

Product Number: **428272**

MDL: MFCD00217740



TEST

Specification

Appearance (Color)	Yellow
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Acid Value	≥ 70.0 meq
Solid Content	88.0 - 92.0 %
Viscosity	700 - 1000 cps

Specification: PRD.0.ZQ5.10000018931

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.