

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

3-Butynyl 2-bromoisobutyrate - 97% (GC)

Product Number:

765104

CAS Number:

1264291-66-8

Formula:


C₈H₁₁BrO₂

Formula Weight:

219.08 g/mol

Storage Temperature:

2 - 8 °C

**TEST****Specification**

Appearance (Color)

Colorless to Faint Brown

Appearance (Form)

Infrared spectrum

Titration by AgNO₃

Purity (GC)

Colorless to Faint Yellow and Colorless to Faint Brown-Yellow and

Liquid

Conforms to Structure

95.5 - 104.5 %

≥ 96.5 %

Specification: PRD.0.ZQ5.10000029639

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.