

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

 β -Glucuronidase from Escherichia coli - aqueous glycerol solution, $\geq 5,000,000$ units/g protein, pH 6.8 (biuret)

Product Number:

G8162

MDL:

MFCD00130629

TEST	Specification
------	---------------

mg protein/ml (Biuret)

1.0 - 10.0

unit/g protein (Biuret)

 ≥ 5000000

One Sigma or modified Fishman unit will liberate 1.0 microgram of phenolphthalein from phenolphthalein glucuronide per hour at pH 6.8 at 37 deg C.

Specification: PRD.1.ZQ5.10000014186

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.