

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 226823
CAS Number: 7783-96-2
MDL: MFCD00003412
Formula: AgI
Formula Weight: 234.77 g/mol

AgI

TEST	Specification
------	---------------

Appearance (Color)	Yellow to Yellow-Green
Appearance (Form)	Powder
Titration by AgNO ₃ % Ag After NaCN Treatment	45.2 - 46.6 %
ICP Major Analysis Confirms Silver Component	Confirmed

Specification: PRD.0.ZQ5.10000030118

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.