

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phorbol 12-myristate 13-acetate - ≥99% (TLC), film or powder

Product Number:

P8139

CAS Number:

16561-29-8

Formula:


C₃₆H₅₆O₈

Formula Weight:

616.83 g/mol

Storage Temperature:

-20 °C

**TEST****Specification**

Solubility (Color)

Colorless to Light Yellow

Solubility (Turbidity)

Clear

10 mg/ml, CHCl₃

Purity (TLC)

≥ 99 %

Recommended Retest Period

2 Years

Specification: PRD.0.ZQ5.10000042226

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.