

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N-(2-Hydroxyethyl)ethylenediamine - 99%

Product Number:**127582**H2NCH2CH2NHCH2CH2OH

CAS Number:

111-41-1

MDL:

MFCD00008170

Formula:

C4H12N2O

Formula Weight:

104.15 g/mol

TEST**Specification**

Appearance (Color)

Colorless

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

Purity (Titration by HCL04)

98.5 - 101.5 %

Purity (GC)

≥ 98.5 %

Specification: PRD.0.ZQ5.10000019058

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.