


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Mellitic acid - 99%

Product Number: M2705
CAS Number: 517-60-2
MDL: MFCD00002469
Formula: C₁₂H₆O₁₂
Formula Weight: 342.17 g/mol

**TEST****Specification**

Appearance (Color)	White to Off-White
Appearance (Form)	Conforms to Requirements
Powder or Crystals	
Infrared spectrum	Conforms to Structure
Titration with NaOH	98.5 - 101.5 %
Anhydrous Basis	
Water (by Karl Fischer)	≤ 5.0 %

Specification: PRD.2.ZQ5.10000005724

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.