

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Olvanil - powder

Product Number: **O0257**
CAS Number: 58493-49-5
MDL: MFCD00673962
Formula: C₂₆H₄₃NO₃
Formula Weight: 417.62 g/mol
Storage Temperature: -20 °C



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
20 mg/ml, DMSO	
Proton NMR spectrum	Conforms to Structure
Carbon	73.2 - 76.2 %
Nitrogen	3.2 - 3.5 %
Purity (TLC)	≥ 98 %

Specification: PRD.1.ZQ5.10000031435

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.