

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

4'-Ethoxyacetophenone – Vetec™ reagent grade

Product Number:

V900621

CAS Number:


1676-63-7

Formula:

C₁₀H₁₂O₂

Formula Weight:

164.20 g/mol



TEST

Specification

Appearance (Colour)

Conforms to Requirements

Appearance (Form)

Conforms to Requirements

Crystals or Crystalline Solid

Infrared spectrum

Conforms to Structure

Purity (GC)

≥ 97.5 %

Specification: PRE.0.ZQ5.10000002515

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.