


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Chloranilic acid

Product Number: **C8136**
CAS Number: 87-88-7
MDL: MFCD00001596
Formula: C₆H₂Cl₂O₄
Formula Weight: 208.98 g/mol



TEST

Specification

Appearance (Color)	Orange to Very Dark Orange
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Color)	Faint Orange to Dark Brown
50 mg/mL, EtOH	
Purity (Titration by NaOH)	≥ 98 %

Specification: PRD.0.ZQ5.10000044333

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.