

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com**Product Specification**

Product Name:

5-Methylthiophene-2-boronic acid pinacol ester 95%

Product Number:

655074

CAS Number:

476004-80-5

MDL:


MFCD05664108

Formula:

C₁₁H₁₇BO₂S

Formula Weight:

224.13 AMU

**TEST****Specification**

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Proton NMR spectra

Conforms to Structure

Purity (GC)

≥94.0 %

Remarks:

Specification Date : 11/22/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.