

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Dexamethasone - powder, γ-irradiated, BioXtra, suitable for cell culture, ≥80% (HPLC)

Product Number: D8893 CAS Number: 50-02-2

MDL: MFCD00064136
Formula: C22H29FO5
Formula Weight: 392.46 g/mol
Storage Temperature: 2 - 8 ℃

HO H<sub>3</sub>C QH QH H<sub>3</sub>C H CH<sub>8</sub>

TEST Specification

Solubility (Color) Colorless

Solubility (Turbidity) Clear

1 mL/VL, EtOH

Sterilized (Gamma Irradiated) Pass

Irradiation dose 25.0 to 35.0 kGy

Purity (HPLC) ≥ 80 %
Cell Culture Test Pass

Specification: PRD.0.ZQ5.10000047195

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.