


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Clobazam

Product Number: C8414
CAS Number: 22316-47-8
MDL: MFCD00079069
Formula: C₁₆H₁₃ClN₂O₂
Formula Weight: 300.74 g/mol



TEST	Specification
Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear
50mg/mL, acetone	
Infrared spectrum	Conforms to Structure
Carbon	62.6 - 65.2 %
Nitrogen	8.7 - 9.7 %
Purity (HPLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000069686

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.