

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com


Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetrabutylammonium bromide - ReagentPlus®, ≥99.0%

Product Number: 193119
CAS Number: 1643-19-2
MDL: MFCD00011633
Formula: C₁₆H₃₆BrN
Formula Weight: 322.37 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Powder or Crystals or Granules
Infrared spectrum	Conforms to Structure
Titration by AgNO ₃	≥ 99.0 %
Solubility (Turbidity)	Clear
50mg/ml, Ethanol	
Solubility (Color)	Colorless to Light Yellow
Recommended Retest Period	-----
20 months	

Specification: PRD.2.ZQ5.10000007186

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.