

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lanosterol - ≥93%, powder

Product Number: L5768 CAS Number: 79-63-0

MDL: MFCD00021108 Formula: C30H50O

Formula Weight: 426.72 g/mol

Storage Temperature: -20 ℃

TEST

Appearance (Color)
Appearance (Form)
Solubility (Color)

Solubility (Turbidity)

50 mg/mL in Chloroform

Infrared spectrum

Carbon Purity (GC)

Specification: PRD.0.ZQ5.10000019162

Specification

White to Off-White

Powder Colorless

Clear

Conforms to Structure 81.5 - 87.4 % > 93 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.