

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Lithium chloride - anhydrous, beads, - 10 mesh, ≥99.9% trace metals basis

Product Number:**449881**

LiCl

CAS Number:

7447-41-8

MDL:

MFCD00011078

Formula:

CLi

Formula Weight:

42.39 g/mol

TEST	Specification
------	---------------

Appearance (Color)

White

Appearance (Form)

Beads

ICP Major Analysis

Confirmed

Confirms Li Component

Trace Metal Analysis

≤ 1000.0 ppm

Purity

Conforms

≥99.9% Based On Trace Metals Analysis

Specification Date : 05/04/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.