

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Glucagon-Like Peptide 1 Fragment 7-37 human - ≥96% (HPLC)

Product Number:

G9416



Formula:

C₁₅₁H₂₂₈N₄₀O₄₇

Formula Weight:

3,355.71 AMU

TEST

Specification

Appearance (Form)

Powder

Appearance (Color)

White

Solubility (Solvent)

Water or 1% Acetic Acid

Solubility (Conc)

1.00 - 1.04 mg/ml

Solubility (Color)

Colorless

Solubility (Turbidity)

Clear

Molecular Weight (MS)

3352.6 - 3358.6 Da

Peptide Content (AAA)

≥ 65.0 %

Purity (HPLC)

≥ 96.0 %

Recommended Retest Period

5 Years

Specification: PRE.0.ZQ5.10000015068

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.