

sigma-aldrich.com

3050 Spruce Street, Saint Louis, MO 63103, USA Website: www.sigmaaldrich.com Email USA: techserv@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Lutetium (III) oxide - 99.99% trace metals basis

Product Number:	203661	
CAS Number:	12032-20-1	
MDL:	MFCD00011100	
Formula:	Lu2O3	
Formula Weight:	397.93 g/mol	

Lu_2O_3

TEST	Specification	
Appearance (Color)	White	
Appearance (Form)	Powder	
Complexometric EDTA	86.6 - 89.3 %	
(% Lu)		
ICP Major Analysis	Confirmed	
Confirms Lutetium Component		
Purity	Meets Requirements	
99.99% Purity Based On Rare Earth Analysis		
Trace Rare Earth Analysis	<pre>< 150.0 ppm</pre>	

Specification: PRD.1.ZQ5.10000023188

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.