

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

SigmaMAb Antibody Drug Conjugate (ADC) Mimic - Antibody Cysteine-Fluorophore Conjugate Standard

Product Number: MSQC8

Storage Temperature: -20 ℃

## TEST Specification

Appearance (Color)

Appearance (Form)
Solubility (Turbidity)

0.5mg/0.1ml H2O

Solubility (Color)

Vial Content

Protein per vial

Mass Spectrometry

Average DAR - Reduced SEC-MS

Minimal Artifactual Modifications

/Adducts - Reduced SEC-MS

Specification: PRD.0.ZQ5.10000080287

White to Off White

Pow der

Clear to Slightly Hazy

Colorless to Light Yellow to Light Brown

> 500 UG

3.2 - 4.8

< 5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.