

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Disialoganglioside GD1b from bovine brain - ~95%, lyophilized powder

Product Number:

G8146

TEST**Specification**

Appearance (Form)	Powder or Lyophilized Solid
Appearance (Colour)	White to Light Yellow
Solubility (Solvent)	CHCl ₃ /MeOH 1/1
Solubility (Conc)	9.80 - 10.20 mg/ml
Solubility (Turbidity)	Clear to Slightly Hazy
Solubility (Color)	Colorless to Faint Yellow
Purity (TLC)	≥ 93.00 %
Recommended Retest Period	-----
3 Years	

Specification: PRE.0.ZQ5.10000001157

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.