

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Calcium Channel (α_1C Subunit) (L-type of Voltage-gated Ca^{2+} Channel) antibody produced in rabbit - affinity isolated antibody

Product Number:

C1603

MDL:

MFCD02095952

TEST**Specification**

Formulation	Conforms
Lyophilized from phosphate buffered saline, pH 7.4, containing 1% BSA, and 0.05% sodium azide	
Protein Content	0.3 - 1.0 mg/ml
after reconstitution	
Working Dilution	Conforms
by Direct Immunoblotting	
1:200 using rat brain and heart membranes, anti-rabbit IgG peroxidase conjugate and detection by ECL	

Specification: PRD.0.ZQ5.10000035273

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.