

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Conalbumin from chicken egg white - Substantially iron-free

Product Number: **C0755**
CAS Number: 1391-06-6
MDL: MFCD00130856

Storage Temperature: -20 °C

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Turbidity) 25 mg/mL, H ₂ O	Clear to Slightly Hazy
Water (by Karl Fischer)	≤ 7 %
Fe by ICP	≤ 0.02 %
Iron	
Agarose Electrophoresis	≥ 98 %

Specification: PRD.1.ZQ5.10000028186

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.