

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Titanium - evaporation slug, diam. x L 6.3 mm x 6.3 mm, ≥99.99% trace metals basis

**Product Number:** 433667  
**CAS Number:** 7440-32-6  
**MDL:** MFCD00011264  
**Formula:** Ti  
**Formula Weight:** 47.87 g/mol

Ti

TEST	Specification
Appearance (Color)	Conforms to Requirements
Grey to Silver	
Appearance (Form)	Conforms to Requirements
Pellets or Slugs	
ICP Major Analysis	Confirmed
Confirms Titanium Component	
Trace Metal Analysis	≤ 100.0 ppm
Purity	Meets Requirements
> = 99.99% Based on Trace Metals Content	
Size	Confirmed
6.3mm Diameter x 6.3mm Length, Approximately	

Specification: PRD.0.ZQ5.10000031359

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.