

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Tumor Necrosis Factor- α antibody produced in goat - affinity isolated antibody, lyophilized powder

Product Number:

T0938

TEST	Specification
Formulation Lyophilized from a 0.2 μ m-filtered solution of PBS with 5% Trehalose	Conforms
Bioactivity ND50 (TNF-Alpha neutralization, L-929 cells)	0.1 - 0.4 μ g/ml
Working Dilution by Capture EI 0.2-0.8 μ g/mL detects recombinant mouse TNF-Alpha	Conforms
Working Dil. by Indirect Immun 0.1 μ g/mL detects 0.5 ng of recombinant mouse TNF-Alpha/lane under non-reducing and reducing conditions	Conforms
Endotoxin Assay < 0.1 EU/ μ g antibody	Conforms

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.