


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Mesityl-2,5,6,7-tetrahydropyrrolo[2,1-c][1,2,4]triazol-4-ium chloride - 97%

Product Number: **688487**
CAS Number: 862893-81-0
Formula: C₁₄H₁₈ClN₃
Formula Weight: 263.77 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Light Yellow and Faint Beige to Beige
Appearance (Form)	Crystal to Powder
Infrared Spectrum	Conforms to Structure
Carbon	61.5 - 66.0 %
Nitrogen	15.3 - 16.5 %

Specification: PRD.0.ZQ5.10000052368

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.