

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
β-Lactoglobulin A from bovine milk - ≥90% (PAGE)

Product Number: L7880
CAS Number: 50863-92-8
MDL: MFCD00130636

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear to Slightly Hazy
10 mg/mL H ₂ O	
Nitrogen	14 - 16 %
Purity (by Electrophoresis)	≥ 90 %
polyacryl-amide electrophoresis	

Specification: PRD.0.ZQ5.10000067971

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.