

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Strontium peroxide - 98%

Product Number: 415200
CAS Number: 1314-18-7
MDL: MFCD00084686
Formula: O₂Sr
Formula Weight: 119.62 g/mol



TEST	Specification
Appearance (Color)	White to Beige
Appearance (Form)	Powder
X-Ray Diffraction	Conforms to Structure
Complexometric Titration	
Typically 72% - 73% Strontium	
Titration with KMnO ₄	
Typically 12% - 13% Available Oxygen	
Specification: Report Result	

Specification: PRD.2.ZQ5.10000068261

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

