

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Platinum - foil, thickness 0.1 mm, 99.9% trace metals basis

**Product Number:** 349380  
**CAS Number:** 7440-06-4  
**MDL:** MFCD00011179  
**Formula:** Pt  
**Formula Weight:** 195.08 g/mol

**Pt**

TEST	Specification
------	---------------

Appearance (Color): Silver	Conforms
Appearance (Form): Foil	Conforms
Thickness: 0.09 - 0.11mm	Conforms
ICP Major Analysis	Confirmed
Confirms Platinum Component	
Purity	Conforms
99.9% Based on Trace Metals Analysis	
Trace Metal Analysis	≤ 1500.0 ppm

Specification: PRD.0.ZQ5.10000039787

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.