

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Petroleum ether – ACS reagent

Product Number: **320447**
CAS Number: 8032-32-4
MDL: MFCD00081849

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Color Test	≤ 10 APHA
Refractive index at 20 °C	1.358 - 1.368
Residue on Evaporation	≤ 0.001 %
Acidity	Pass
Initial Boiling Point	≥ 35 °C
Final Boiling Point	≤ 60 °C
Meets ACS Requirements	Current ACS Specification

Specification: PRD.2.ZQ5.10000036441

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.