

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
N-TER™ Nanoparticle siRNA Transfection System

Product Number: **N2913**

Storage Temperature: -20 °C

TEST	Specification
Performance	Conforms
NLT 70% knockdown and NLT 70% viability using 30nM or less GAPDH siRNA in Hela cells	
Sterility	Negative

Specification: PRD.1.ZQ5.10000095603

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.