

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
CompoZr® Targeted Integration Kit – AAVS1 – Transgene Integration in the Human Genome

Product Number: CTI1

Storage Temperature: -20 °C

TEST	Specification
Suitability	Suitable
Suitable for targeted integration.	
Concentration	0.4 - 0.8 µg/µl
Component M7947	
Concentration	1.5 - 2.0 µg/µl
Component Z3027	
DNA Concentration	20 - 50 ng/µl
Component I2410	

Specification: PRD.3.ZQ5.10000005461

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.