

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium sulfite – anhydrous, Redi-Dri™, reagent grade, ≥98%

Product Number: 901916
CAS Number: 7757-83-7

Formula: Na₂SO₃
Formula Weight: 126.04 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Crystals
Sodium	Positive
Sodium ID by ICP	
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
100 mg/mL, H ₂ O	
Titration with Na ₂ S ₂ O ₃	≥ 98 %
Recommended Retest Period	-----
2 Years	

Specification: PRD.5.ZQ5.10000088528

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.