

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Coumermycin A1

Product Number: **C9270**
CAS Number: 4434-05-3
MDL: MFCD00057312
Formula: C₅₅H₅₉N₅O₂₀
Formula Weight: 1,110.08 g/mol



TEST	Specification
Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
50 mg/ml DMSO	
Infrared spectrum	Conforms to Structure
Purity (TLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000059016

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.